

OFFICE OF THE MISSION DIRECTOR NATIONAL HEALTH MISSION, ASSAM

SAIKIA COMMERCIAL COMPLEX, CHRISTAN BASTI, G.S ROAD, GUWAHATI-781005 PH.NO: 0361-2363062; TELE FAX: 0361-2363058

Website: https://nhm.assam.gov.in

e-mail: misnrhm.assam@gmail.com

TENDER NO: NHM/18019/1/2019-PROC-NHM/ 6 7 29

DATE: 21/6/2119

CORRIGENDUM NO.1 NOTICE INVITING TENDER (NATIONAL COMPETITIVE BIDDING)

This is with reference to the TENDER NO: NHM/18019/1/2019-PROC-NHM/6153, DATE: 18/6/2019.


The following amendments in the tender should be taken note of prior to submission of bids.

1. Clause 1.2.1 Earnest Money Deposit (EMD) should be read as bellow EMD shall be @ Rs 1,20,000.00 in the form of Demand Draft/Bankers Cheque in favour of "STATE HEALTH SOCIETY, ASSAM".

Without minimum EMD, a bid will not be considered at all.

2. ANNEXURE - VII-B <u>DECLARATION ON MANUFACTURING CAPACITY & SPARABLE CAPACITY:</u> should be read as <u>"MODIFIED ANNEXURE-VII-B(Rev)</u>

SI. No	Name of the item with specification	Unit	Estimated Annual Requirement	Annual Production Capacity	Annual Separable Capacity for Assam
1	Dengue NS 1 Antigen Detection Kit by Elisa	1 Unit = A single kit of 1 X 96 test	109 Kit		
2	DELTAMETHRINE 2.5% FLOW	Bottle of 1 Liter each	12900 Bottle		
3	Malathion	Drum of 250 kg	64 Drum (16000 kg)	,	



298

3. ANNEXURE— XIII: ITEMS WITH QUANTITIES: should be read as "MODIFIED ANNEXURE— XIII (Rev): ITEMS WITH QUANTITIES" and the list of item with quantity is as follows:-

- 84

LIST OF ITEMS WITH QUANTITIES

Tender Ref.	No:	Date:
i ciidei iteii		

SI. No	Name of item	Dosage form and strength	Pack Size	Tentative Quantity required for 2 year	Price Quoted/ Not Quoted
1	Dengue NS 1 Antigen Detection Kit by Elisa	 The Elisa Kit should be designed for qualitative detection on Dengue NS1 antigen of all 4 Dengue serotypes in human serum. The kit should be provided with the following materials and reagents: Anti –NS1 Antibody Coated Breakway Microwells (12 X 8 =96). Desiccant should be provided for storing the unused microwells which are to be resealed immediately. Horseradish peroxides conjugated Anti-NS1 monoclonal antibody with preservative. Chromogenic substrate in buffer. Positive control in the form of recombinant antigen with preservatives and antibodies. Negative control in the form of confirmed negative human serum with preservatives and antibodies. Calibrators in the form of recombinant antigen with preservatives and antibodies. Calibrators in the form of secombinant antigen with preservatives and antibiotics. Sample diluents. Wash Buffer. Wash Buffer. Wash Buffer. Wash Buffer. Mash Buffer. Positive control in the form of recombinant antigen with preservatives and antibiotics. Mash Buffer. Mash Buffer. Mash Buffer. Positive control in the form of recombinant antigen with preservatives and antibiotics. Mash Buffer. Positive control in the form of recombinant antigen with preservatives and antibiotics. Positive control in the form of recombinant antigen with preservatives and antibiotics. Positive control in the form of recombinant antigen with preservatives and antibiotics.	1 Unit = A single kit of 1 X 96 test	218 Kit	







		test for deter antigen by I between 2- 4) The Elisa kit NS1 should 90% and a strangen Br. PCR as 5 5) The kit should least 6 mon ambient terms.	quired for performing the ection of Dengue NS1 Elisa should range 4 hours. It for detection of Dengue have a sensitivity of > specificity of >95% taking the gold standard. Id have a self – life of at ths when stored at an inperature of 2°C -8°C hould be under cold		
		<u>IS No</u> 6940:1982	Title Methods of test for pesticides and their formulations (first revision)		
	DELTAMETHR	8190:1988(Part)	Requirement of packing of pesticides: Part 2 liquid pesticides (Second edition)		
2	INE 2.5% FLOW: the product is intended for use in impregnation of Bed Nets(ITBN)	10627:1983 12005:1987	Methods for sampling of pesticide formulations Deltamethrine Technical	1 Bottle =1 li	12900 bottle
		Deltameth IS 12005: solvent(s) antifoami agents and b. Descripti The mater white to compare the solution of th	rial shall consist of nrine Technical confirming to 1987 dissolved in suitable to together with thickening ng, wetting and dispersing d with suitable stabilizers. on rial shall be in the form of cream colored liquid. On with water, it shall be suitable		







c. Deltamethrine Content

The nominal Deltamethrine content shall be declared, when tested by the method prescribed in Appendix A of IS 12005:1987, the observed Deltamethrine content percent by mass, (m/m), shall not differ from declared nominal value by more than the tolerance limits given below.

Nominal value percent	Tolerance
limit	
Up to 9	+10
	-5
Above 9 and below 50	+-5
50 and above	+5
	-3

d. Dispensability ion water
The material shall pass the test
prescribed as below:Fill up a 100ml measuring cylinder with
95ml semi hard water and add 5ml of the
sample. Stopper the cylinder and shake
thoroughly thirty times. Allow the
cylinder to stand at 27 Degree C for one
hour.

The sample is considered to be passing the test if it remains dispersed in water and any separation at the top or sedimentation at the bottom is not more than 2ml.

e. Acidity

When tested by the method prescribed in 13.5 of IS6940:1982.

Acidity (as H2SO4) shall not be more than 0.5 percent by mass (m/m).

iv) PACKING

The material shall be packed in one liter HDPE containers which shall be encased in fibre board boxes or wooden boxes. It shall conform to the general requirements given in IS 8190 (Part 2): 1988

v) MARKING

The containers shall bear legibly and inedible the following information:

Name of the material:

Name of the manufacturer:

Batch Number:

Date of Manufacturer:

Date of Expiry:

Nominal Deltamathrine Count (m/m) as registered under the Provisions of Insecticides Act, 1968 and the rules framed thereunder; Volume in ml or L; and Any other information required under the







	Insecticides Act,1968 and the Rules framed thereunder and Standards of Weights and Measures (Packaged Commodities) Rules, 1977 SAMPLING When bulk material is offered for inspection, representatives samples of the material shall be drawn as prescribed in IS 10627: 1983, and tested within 90 days of its manufacturer. The criteria of conformity shall be as given in IS 10627:1983. When the material is offered for inspection after 90 days of its manufacturer, sampling shall be done as prescribed in IS 10627-1983, however, the criteria for conformity of the material shall be the limits of tolerances given under the heading "Deltamathrine Content" of this specification. SHELF LIFE: - Shelf life from manufacturing day to expiry date should be at least 2 years and it should not pass more than 1/6th of their effective life from the date at the time, the material is supplied.		
Malathion	Description: The material shall be in the form of colorless to amber liquid free from extraneous impurities, sediments or added modifying agents. Purity (% by mass Min.) 95 Acidity as H2SO4 (% by mass Max) 0.40 Moisture (% by mass Max) 0.10 Material insoluble in acetone (% by mass Max) 0.50 Specify Gravity at 25'C / 25' Min 1.23 Material Conforms to standard IS: 1832-1978	64 Drum (16000 kg)	

a

<u>NOTE:</u> The quantity against each item is tentative and subject to variations. The quantity may increase or decrease at the time of placing order.







All other terms & Conditions of the above referred tender shall remain unchanged.

P

Mission Director, NHM, Assam

Memo NO: NHM/18019/1/2019-PROC-NHM/ 6 + 30 - 35 Copies to the:

DATE: 21/6/19

- 1. Principal Secretary to the Govt. of Assam, Health & F.W. Deptt. for kind information.
- 2. Director of Health Services, Hengerabari, Guwahati -for information.
- 3. Director Medical Education, Six Mile, Khanapara, Guwahati -for information.
- 4. MIS Manager, NHM, with a request to upload the Corrigendum Notice in the NHM Website.
- 5. SME, NHM, with a request to publish the Corrigendum Notice in leading national & local newspapers.
- 6. P.S to the Honb'le Minister, H &FW, Govt. of Assam for kind appraisal of Honb'le Minister.

Executive Director, NHM, Assam